



2236 3M RED DOT TRACE PREP 18 MM X 5 M 10.7PO X 196 PACK OF 10

2236 3M Red Dot Trace Prep 18 mm x 5 m 10.7po x 196 Pack of 10

Highlights


Red Dot™ Trace Prep abrades skin prior to electrode placement

Gently removes top non-conductive layer of skin

Improves trace quality

3/4" wide, 196" long

Convenient case of 36 rolls



SKU: N/A

Price:

Stock: N/A

Category: [PPE](#)

1

PRODUCT DESCRIPTION

2236 3M Red Dot Trace Prep 18 mm x 5 m 10.7po x 196 Pack of 10

Highlights

Red Dot™ Trace Prep abrades skin prior to electrode placement

Gently removes top non-conductive layer of skin

Improves trace quality

3/4" wide, 196" long

Convenient case of 36 rolls

