



3F020120 TRAVELER® L4 MANUAL FOLDING WHEELCHAIRS

3F020120 Traveler® L4 Manual Folding Wheelchairs

Available in 16", 18", 20" seat widths. Also 18" and 20" widths with 18" seat depth

Maximum Weight Capacity 250 lb (300 lb on models with 18" seat depth) - EVENLY DISTRIBUTED

Five (5) year limited warranty



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DURABLE MEDICAL](#), [Wheel Chairs](#)

Tags: [traveler wheelchairs](#), [wheelchairs](#)

1

PRODUCT DESCRIPTION

3F020120 Traveler® L4 Manual Folding Wheelchairs

The Traveler L4 is a high strength, lightweight Medicare approved K0004 wheelchair designed for both the home and rehab markets. Weighing less than 34 lbs, it features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers flip-back desk and full length arm options for easy lateral movement. Precision bearings and non-marking polyurethane rear tires also expand the longevity of the wheelchair.

Standard Features

- Available in 16", 18", 20" seat widths. Also 18" and 20" widths with 18" seat depth
- Quick release axle is available for some models
- Weighs less than 34 lb without front rigging



- Hammertone paint provides a long lasting and durable finish
- Carbon steel frame provides strength and durability
- Black nylon upholstery with standard chart pocket on back
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Flip-back arms
- Height adjustable flip-back arms on selected models
- Height adjustable back
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb (300 lb on models with 18" seat depth) - EVENLY DISTRIBUTED
- Meets Medicare code K0004
- This product is not made with natural rubber latex

Warranty

Five (5) year limited warranty