




4299C PATIENT PRIVACY SCREEN 3 PANELS

4299C Patient Privacy Screen 3 Panels

Each Panel: 27.5" x 56" (700 mm x 1,420 mm)

Overall Height: 70.8" (1,800 mm)

Weight: 6.75 kg



SKU: N/A

Price:

Stock: Out of stock

Categories: [PATIENT CARE](#), [Privacy Screens](#)

Tags: [patient screen](#), [privacy screen](#)

1

PRODUCT DESCRIPTION

4299C Patient Privacy Screen 3 Panels

4299C (4) Replacement Casters and (4) Tips

6mil white vinyl panels are flame-retardant*

Easy to set up and to dismantle for storage

Vinyl panel is washable and replaceable on snap-out rods

Frame is 1" (25mm) anodized aluminum tubing

2" (51mm) spherical casters allow for easy movement

To clean vinyl panels, use mild liquid soap and water, rinse.

*NFPA 701 Test Method 1-2010 Edition

Specifications



Each Panel: 27.5" x 56" (700 mm x 1,420 mm)

Overall Height: 70.8" (1,800 mm)

Weight: 6.75 kg