



01690-200 WELCH ALLYN SURE TEMP ELECTRONIC THERMOMETER (MODEL 690)

01690-200 Welch Allyn Sure Temp Electronic Thermometer (MODEL 690)

Item No: 01690-200

Large LCD displays temperatures in Fahrenheit or Celsius

Last temperature recall button



SKU: N/A

Price:

Stock: Available on backorder

Categories: [DIAGNOSTIC INSTRUMENTS](#), [Thermometry](#)

Tag: [thermometer](#)

1

PRODUCT DESCRIPTION

01690-200 Welch Allyn Sure Temp Electronic Thermometer (MODEL 690)

Item No: 01690-200

Help your facility improve staff productivity and patient safety with the SureTemp Plus 690 Electronic Thermometer. The SureTemp 690 Thermometer is fast-capturing approximately 4- to 6-second oral, 10- to 13-second rectal or pediatric axillary and 12- to 15-second adult axillary temperatures. Featuring the technology that sets the standard for accuracy, you can have confidence the SureTemp 690 will deliver accurate, repeatable readings on patients of all age groups. The removable/interchangeable oral/axillary and rectal probes and probe wells are designed to help you reduce the risk of cross contamination, for staff satisfaction and patient safety.



Features

Fast: quickly captures accurate readings on all patient types

Large LCD displays temperatures in Fahrenheit or Celsius

Convenient storage for 25 probe covers

Last temperature recall button

Optional security wall mount

Takes approximately 6,000 readings on 3 AA batteries

Interchangeable, removable wells to help reduce the risk of cross-contamination