




058-701 BIODEX TABLE ECHOCARDIOGRAPHY ECONO 115V

058-701 Biodex Table Echocardiography

Econo 115V

Item No: 058-701

The Econo Echocardiography table is based on the same platform as the Echo Pro.



SKU: N/A

Price:

Stock: Available on backorder

Categories: [IMAGING TABLES](#), [Ultrasound Table](#)

Tags: [imaging table](#), [ultrasound table](#)

1

PRODUCT DESCRIPTION

058-701 Biodex Table Echocardiography

Econo 115V

Item No: 058-701

The Econo Echocardiography table is based on the same platform as the Echo Pro. It has the same ergonomic features, keeping the safety of the sonographer in mind. To make it a more cost-effective alternative, it is offered without Trendelenburg positioning or the central floor-locking system.

The barrier between sonographer and patient is eliminated with the Econo Echocardiography Table, featuring an ergonomic design, which helps prevent musculoskeletal injury.



With the capability to release and return from either side of the table, the cardiac cushion drops down to reveal the industrys largest scanning cutout providing open access to the patients left thorax area for an unobstructed, apical approach.

The Econo Echocardiography Table has many optional items to further improve sonographer comfort and safety. Choose what best suit your departments specific needs. Side rails for closer patient access, an articulating scanning arm board, adjustable from 0 to 130-degrees and a headrest for easier thyroid and carotid access are all available.

Suitable Procedures:

Abdomen

Aorta

Thyroid

Breast

Carotid

Upper Extremity Arterial*

Renal Transplant

SMA

Stress Echo

Renal

Pelvis

Testicular

Musculoskeletal

Upper Extremity Venous*

Renal Artery

Abdominal Doppler

Echocardiography

Transesophageal Echo

*Requires optional vascular arm board.

All motorized tables are ETL product listed to the UL 60601-1 Medical Electrical Equipment Safety Standard.