



## 12022110 SPECIALTY WHEELCHAIRS

12022110 Specialty Wheelchairs

W/C RSC 14 1/4" LoBack with 24" Rear wheel 250 lb E&J

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care.



**SKU:** N/A

**Price:**

**Stock:** 5 in stock (can be backordered)

**Categories:** [BATHROOM SAFETY](#), [Commodes](#)

**Tags:** [rehab commode](#), [shower commode](#)

1

## PRODUCT DESCRIPTION

12022110 Specialty Wheelchairs

12022110 W/C RSC 17 1/4" HiBack with 24" Rear wheel 250 lb E&J

The versatile E&J Folding Rehab Shower Commode chair has earned the reputation of being the premier self-propelled chair for personal hygiene care. This chair features an overall narrow width and a corrosion-resistant stainless steel frame. The chair can be positioned over a standard toilet or used as a portable, self-contained commode.

Standard Features

- Completely removable full length arms for ease of transfer
- Swingaway, detachable footrests permit close approaches to beds, toilets and wheelchairs
- Recessed seat and frame side contours allow easy access for toileting activities while seated



- Available with 24" composite wheels with integral handrims and solid rubber tires or 5" rear casters on models with 14 1/4" & 17 1/4" back height
- Solid nylon tires on models with 16" & 19" back height adding more durability
- Removable, hygienic plastic pail and lid are included
- Cushioned vinyl seat is sealed against moisture penetration and vented for optimal foam compression
- Non-corrosive composite footplates
- Maximum weight capacity 250 lb EVENLY DISTRIBUTED on models with 14 1/4" & 17 1/4" back height
- Maximum weight capacity 300 lb EVENLY DISTRIBUTED on models with 16" & 19" back height
- Chair comes standard with oval opening seat C-shaped seat is an optional item
- Meets Medicare code E0165

#### Warranty

Two (2) year limited warranty