



1280 GOLDNAMER EAR BASIN

1280 Goldnamer Ear Basin

conforms to any ear size. One-piece, seamless white plastic. Lightweight. Specially molded to fit under the patient's ear. Autoclavable up to 250° F and/or cold sterilizable.



SKU: N/A

Price:

Stock: 3 in stock

Category: [PPE](#)

Tags: [Ear Basin](#), [Goldnamer Ear Basin](#)

1



PRODUCT DESCRIPTION

1280 Goldnamer Ear Basin

conforms to any ear size. One-piece, seamless white plastic. Lightweight. Specially molded to fit under the patient's ear. Autoclavable up to 250° F and/or cold sterilizable.

1