



140 IV TRAINING ARM AND HAND SET

140 IV Training Arm and Hand Set

These training aids were designed for practice in giving IVs. They permit detailed exercises in venipuncture so that the technique can be mastered on a realistic model before actually administering to patients.



SKU: N/A

Price:

Stock: Available on backorder

Category: [PPE](#)



ADDITIONAL INFORMATION

IV Training Arm and Hand Set Light I Dark

Color [Dark](#), [Light](#)

PRODUCT DESCRIPTION



140 IV Training Arm and Hand Set

These training aids were designed for practice in giving IVs. They permit detailed exercises in venipuncture so that the technique can be mastered on a realistic model before actually administering to patients. The students can palpate the arm and hand, which are made of supple, resilient plastic with a human feel. The vein is located, the needle introduced, and blood is withdrawn or fluid injected. The models can be used repeatedly. Comes complete with blood powder and carrying trays. No needle included. Three-year warranty.

[please check out our other PPE products.](#)

ADDITIONAL INFORMATION

IV Training Arm and Hand Set

Light I Dark

Color

[Dark](#), [Light](#)