




2188 GRAVITY INCLINOMETER

2188 Gravity Inclinometer
INCLINOMETER GRAVITY



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [gravity inclinometer](#), [inclinometer](#)

1

PRODUCT DESCRIPTION

2188 Gravity Inclinometer
INCLINOMETER GRAVITY

Measures range of motion of any joint. Moveable dial for manual "O" setting. Calibrated in 1" increments

1