



22" SPRAGUE RAPPAPORT-TYPE STETHOSCOPE

22" Sprague Rappaport-Type Stethoscope

602 BLK

602BL BLU

602R RED

602TL TEAL

602LB LT BL

602MAG MAGENTA

602L LAV

602BUR BURG

600ID Stethoscope ID Tag

DP600SM Diaphragm, Small (Pediatric), for Sprague-Rappaport Type

SR610 Accessory Kit



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DIAGNOSTIC INSTRUMENTS](#), [Stethoscopes](#)

Tags: [rappaport-type stethoscope](#), [sprague stethoscope](#), [stethoscope](#)

ADDITIONAL INFORMATION

Color 107, 251, 105, 250, 249, 85, 242

Item No. 4 1691, 1690



VARIATIONS

Image	SKU	Price	Description	Color	Item No. 4
				Black	SR610
				Black	600ID
				Red	Any Item No. 4
				Burgundy	Any Item No. 4



Image	SKU	Price	Description	Color	Item No. 4
				Black	Any Item No. 4
				Slate Blue	Any Item No. 4
				Magenta	Any Item No. 4
				Lavender	Any Item No. 4



Image	SKU	Price	Description	Color	Item No. 4
				Blue	Any Item No. 4

PRODUCT DESCRIPTION

22" Sprague Rappaport-Type Stethoscope

602 STETH SPRAGUE-RAP TYPE 22 BLK

602BL STETH SPRAGUE-RAP TYPE 22 BLU

602R STETH SPRAGUE-RAP 22" RED

602TL STETH SPRAGUE-RAP 22" TEAL

602LB STETH SPRAGUE-RAP 22" LT BL

602MAG STETH SPRAGUE-RAP 22" MAGENTA

602L STETH SPRAGUE-RAP 22" LAV

602BUR STETH SPRAGUE-RAP 22" BURG

600ID Stethoscope ID Tag

DP600SM Diaphragm, Small (Pediatric), for Sprague-Rappaport Type

SR610 Accessory Kit

Item 602 Series

* Five interchangeable chestpieces (Two diaphragms, Three bells)
facilitates detection of high, medium, and low frequency sounds

* Ideal for infant, pediatric and adult use

* Includes three ear tip sets to maximize comfort

* Steel binaural

* Heavy-walled 22" tubing decreases extraneous sounds

* Accessory pouch included

* Latex free



ADDITIONAL INFORMATION

Color	107, 251, 105, 250, 249, 85, 242
Item No. 4	1691, 1690