



240 HOME BLOOD PRESSURE KIT WITH SEPARATE STETHOSCOPE

240 Home Blood Pressure Kit with Separate Stethoscope

240 Adult

240C Child

240X Large Adult



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Blood Pressure Sphygmomanometer](#),
[DIAGNOSTIC INSTRUMENTS](#)

Tags: [blood pressure](#), [blood pressure kit](#)

ADDITIONAL INFORMATION

size age	144, 145, 261
-----------------	---------------

VARIATIONS



Image	SKU	Price	Description	size age
				adult
				child
				Large Adult

PRODUCT DESCRIPTION

240 Home Blood Pressure Kit with Separate Stethoscope

240 Adult

240C Child

240X Large Adult

- Cotton cuff and latex-free two (2) tube bladder
- Self-adjusting D-Bar cuff
- Pressure measurement from 20 mmHg to 300 mmHg
- Ten (10) year calibration warranty



- Zippered carry bag is included

Cuff Class Limb Range (inches) Bladder (inches) Cuff Dimensions (inches)*

Newborn 3.75" - 5.25" 1" x 3.55" 7" W x 2" H

Infant 5" - 7.5" 2.35" x 4.72" 10.5" W x 3" H

Child 7" - 10.5" 3.55" x 5.90" 13.5" W x 4" H

Adult 9.5" - 15" 4.72" x 8.66" 20.5" W x 5.5" H

LG Adult 13" - 19" 5.71" x 12.41" 25" W x 6.5" H

Thigh 17" - 26" 7.09" x 14.57" 29.5" W x 8.5" H

*: These are subjective measurements and can vary $\pm 0.25"$

ADDITIONAL INFORMATION

size	age	144, 145, 261
-------------	------------	---------------