



3227 BED PANS

3227 Bed Pans

Adult: 14" x 11 3/8".

Weight capacity: 500lbs



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [MEDICAL, SURGICAL and PHYSICIAN, Stainless Steelware](#)

Tags: [bed pans, patient pans](#)

1

PRODUCT DESCRIPTION

3227 Bed Pans

Adult: 14" x 11 3/8".

Front: 3 1/4".

Back: 1 15/16"

Weight capacity: 500lbs

Seamless, sanitary, highly-polished 18-8 stainless steel ware

Reusable

Autoclavable

1