




3233 DRESSING JARS

3233 Dressing Jars
6 7/8" x 4 7/8". Capacity: 2 1/8 qt.



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#),
[Stainless Steelware](#)

Tag: [dressing jars](#)

1

PRODUCT DESCRIPTION

3233 Dressing Jars
6 7/8" x 4 7/8". Capacity: 2 1/8 qt.
18-8 (304) Stainless Steel with seamless construction.

1