



3397 DENTURE PLATE BRUSH

3397 Denture Plate Brush
3397 BRUSH DENTURE PLATE
24 / BX



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Personal Care/Daily Living Aids](#), [Sundries](#)

Tags: [denture brush](#), [denture plate brush](#)

1

PRODUCT DESCRIPTION

3397 Denture Plate Brush
3397 BRUSH DENTURE PLATE
24 / BX

Polypropylene white handle 3 5/8" long

Bristled on 2 sides with nylon filaments

One side with 4 rows - 16 tufts with contour trimming for easy reach into crevices in the anatomy of the denture

Other side with 4 rows - 22 tufts with serrated trimming

Can withstand temperature changes from 0° - 180°F

Individually packaged

