



## 3831 VINYL EMERGENCY CADAVER BAG

3831 Vinyl Emergency Cadaver Bag  
3831 36" x 90"



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [infection control](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

**Tags:** [cadaver bag](#), [infection control](#)

1

## PRODUCT DESCRIPTION

3831 Vinyl Emergency Cadaver Bag  
3831 36" x 90"

Six gauge PVC pouch, double-sealed at both ends for extra protection. Full length, center zipper closure. Water resistant.

Thickness - 6MIL (.15MM)

1