




## 3E010100 TRAVELER® SE TRAVELER SE, 18X16 FIXED FULL ARM, SWINGAWAY FOOTREST

3E010100 Traveler SE Traveler SE, 18x16 Fixed Full Arm, Swingaway Footrest

The Traveler SE incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction. This wheelchair comes in a hammertone finish and provides the most effective combination K0001/K0002 wheelchair available.



**SKU:** N/A

**Price:**

**Stock:** 4 in stock

**Categories:** [DURABLE MEDICAL](#), [Wheel Chairs](#)

**Tag:** [wheelchairs](#)

1

### PRODUCT DESCRIPTION

3E010100 Traveler SE, 18x16 Fixed Full Arm, Swingaway Footrest

The Traveler SE incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction. This wheelchair comes in a hammertone finish and provides the most effective combination K0001/K0002 wheelchair available.

Standard Features:

- Available in 16", 18" and 20" seat widths, with 16" seat depth and 20" seat width with 18" seat depth
- Black vinyl upholstery with standard chart pocket on back
- Industry compatible tool-free front riggings with heel loops standard on swingaway footrests
- Hammertone paint provides a long lasting and durable finish



- Two-position axle and caster fork provide 2" seat height adjustment for a hemi position
- Maximum Weight Capacity 250 lb -EVENLY DISTRIBUTED on 16"x16", 18"x16" and 20"x16" models
- Maximum Weight Capacity 300 lb -EVENLY DISTRIBUTED on 20"x18" models
- Meets Medicare code K0001 and K0002 and K0006 on 20"x18" models
- 16" seat depth models are available for Custom Color Program, restrictions apply

**Warranty:**

One (1) year limited warrantyThe Traveler SE incorporates the legendary Everest & Jennings® wheelchair frame with saddle-joint tubular construction. This wheelchair comes in a hammertone finish and provides the most effective combination K0001/K0002 wheelchair available.