



3F012120 TRAVELER® L3 PLUS MANUAL FOLDING WHEELCHAIRS

3F012120 Traveler® L3 Plus Manual Folding Wheelchairs

Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary

Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED

Five (5) year limited warranty



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DURABLE MEDICAL](#), [Wheel Chairs](#)

Tags: [traveler wheelchairs](#), [wheelchairs](#)

1

PRODUCT DESCRIPTION

3F012120 Traveler® L3 Plus Manual Folding Wheelchairs

Everest & Jennings introduces the Traveler® L3 Plus as the lightweight wheelchair addition to the Traveler® line, a brand respected for performance and durability. The Traveler® L3 Plus lightweight frame weighs less than 36 pounds making it the ideal K0003 Medicare coded wheelchair built for short and long term use as well as rental needs.

This lightweight wheelchair comes standard with non-marking polyurethane tires, industry compatible front rigging and adjustable anti-tippers along with a detachable 2" seat extension with hardware. It features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers flip-back desk, adjustable height and full arm options for



easy lateral movement. The Traveler® L3 Plus is manufactured with the same commitment to quality Everest & Jennings is known for.

Standard Features

- Weighs less than 36 lb without front rigging
- Available in 16", 18", 20" seat widths including 2" removable seat extension kit converting the chair from 16" depth to 18" depth if necessary
- Two position axle and caster fork provides 2" seat height adjustment
- Hammertone finish and carbon steel frame provide strength and durability
- Black nylon upholstery with standard chart pocket on back
- Includes fastener tab under seat connecting both seat and back upholstery to eliminate gap
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Flip-back padded arms available in fixed desk length, adjustable height desk and detachable flip-back full length
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0003

Warranty

Five (5) year limited warranty

*This product is not made with natural rubber latex

<http://catalog.bqmsupplies.com/product/3f012120-traveler-l3-plus-lightweight-wheelchairs/>