



## 3F020120 TRAVELER® L4 LIGHTWEIGHT WHEELCHAIRS

3F020120 Traveler® L4 Wheelchairs

Available in 16", 18", 20" seat widths. Also 18" and 20" widths with 18" seat depth  
Maximum Weight Capacity 250 lb (300 lb on models with 18" seat depth) - EVENLY DISTRIBUTED  
Five (5) year limited warranty



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [DURABLE MEDICAL](#), [Wheel Chairs](#)

**Tags:** [traveler wheelchairs](#), [wheelchairs](#)

1

## PRODUCT DESCRIPTION

3F020120 Traveler® L4 Lightweight Wheelchairs

The Traveler L4 is a high strength, lightweight Medicare approved K0004 wheelchair designed for both the home and rehab markets. Weighing less than 34 lbs, it features a dual position axle, multi-position caster forks for 18 or 20 inch seat height adjustment and offers flip-back desk and full length arm options for easy lateral movement. Precision bearings and non-marking polyurethane rear tires also expand the longevity of the wheelchair.

Standard Features

- Available in 16", 18", 20" seat widths. Also 18" and 20" widths with 18" seat depth
- Quick release axle is available for some models
- Weighs less than 34 lb without front rigging



- Hammertone paint provides a long lasting and durable finish
- Carbon steel frame provides strength and durability
- Black nylon upholstery with standard chart pocket on back
- Anti-tippers come standard
- Tool-free adjustable front rigging
- Non-marking polyurethane tires
- Industry compatible front rigging with heel loops standard on swingaway footrests
- Flip-back arms
- Height adjustable flip-back arms on selected models
- Height adjustable back
- Two position axle and caster fork provides 2" seat height adjustment
- Maximum Weight Capacity 250 lb (300 lb on models with 18" seat depth) - EVENLY DISTRIBUTED
- Meets Medicare code K0004
- This product is not made with natural rubber latex

#### Warranty

Five (5) year limited warranty