



3K010230 ADVANTAGE RECLINER WHEELCHAIRS


3K010230 Recliner Wheelchairs

Available in 16", 18" and 20" seat widths, with 17" seat depth and in 22" seat width with 19" seat depth

Maximum weight capacity 300 lb - EVENLY DISTRIBUTED for 16,18 and 20" models

Maximum weight capacity 450 lb - EVENLY DISTRIBUTED for 22" models

Five (5) year limited warranty



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DURABLE MEDICAL](#), [Wheel Chairs](#)

Tags: [recliner wheelchairs](#), [wheelchairs](#)

1

PRODUCT DESCRIPTION

3K010230 Advantage Recliner Wheelchairs

The Advantage Recliner comes with a higher back and removable full support headrest. The headrest includes an additional neck support that adjusts easily to fit every need. The durable hammertone frame is designed to support up to 300 lb and in the bariatric version up to 450 lb users.

Standard Features

- Available in 16", 18" and 20" seat widths, with 17" seat depth and in 22" seat width with 19" seat depth



- Features manual full reclining back range from 90-170 degrees
- Removable full support headrest
- Detachable padded desk or full arms
- Hammertone finish provides durability and style
- Embossed black leatherette type vinyl upholstery for easy cleaning
- Anti-tippers standard
- Comes with spreader metallic back bar which works as anti-fold device
- Industry compatible tool-free adjustable elevating front rigging
- Standard chart pocket on back
- Maximum weight capacity 300 lb - EVENLY DISTRIBUTED for 16,18 and 20" models
- Maximum weight capacity 450 lb - EVENLY DISTRIBUTED for 22" models
- 16", 18", 20" models meet Medicare code K0001/ E1226
(Recliner Back) + E0966 (Head Rest Extension) + E0971
(Anti-tipping device) + K0195(Elevating Legrests)
- 22" models meet Medicare code K0007/E1226
(Recliner Back + E0966 (Head Rest Extension) + E0971
(Anti-tipping device) + K0195 (Elevating Legrests) + E2201

Warranty

Five (5) year limited warranty