



540 TEACHING STETHOSCOPE

540 Teaching Stethoscope

540 DUAL HEAD Black

Adjustable chrome-plated external spring binaural.



SKU: N/A

Price:

Stock: 10 in stock

Categories: [DIAGNOSTIC INSTRUMENTS](#), [Stethoscopes](#)

Tags: [stethoscope](#), [teaching stethoscope](#)

1

PRODUCT DESCRIPTION

540 Teaching Stethoscope

540 STETH TEACHING DUAL HEAD BLK Black

Excellent acoustical clarity. 22" latex-free tubing with connector allows both student and instructor to listen simultaneously. Anodized aluminum combination dual head chestpiece and chrome-plated brass rotating stem. Adjustable chrome-plated external spring binaural.

1