




5PX10620 PARAMOUNT XD BARIATRIC/HEAVY DUTY WHEELCHAIRS

5PX10620 Paramount XD Wheelchairs

Available in 26" and 30" seat widths with 19 ½" seat depth

Maximum weight capacity 650 lb - EVENLY DISTRIBUTED

Five (5) year limited warranty



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DURABLE MEDICAL](#), [Wheel Chairs](#)

Tags: [paramount XD wheelchairs](#), [wheelchairs](#)

1

PRODUCT DESCRIPTION

5PX10620 Paramount XD Bariatric/Heavy Duty Wheelchairs

The Paramount XD is an alternative for the bariatric user who prefers a manual wheelchair. It features 26" and 30" seat widths as well as a weight capacity of 650 lb. Meets Medicare specifications for K0007 for wheelchairs.

Standard Features

- Available in 26" and 30" seat widths with 19 ½" seat depth
- Black reinforced nylon upholstery
- Standard chart pocket on back upholstery for models after serial number AIS338WC0046-0055.
- Reinforced welds and gussets on frame with double crossbrace



- Two-position axle and caster fork provides 2" seat height adjustment
- 8"x 2" non-marking solid casters and 24" solid rear wheels with chromed handrims
- Detachable padded desk arms for easy transfer
- Adjustable swingaway footrest or elevating legrest with extra wide metallic footplate
- Maximum weight capacity 650 lb - EVENLY DISTRIBUTED
- Meets Medicare code K0007

Warranty

Five (5) year limited warranty