



602BK DELUXE SPRAGUE-RAPPAPORT TYPE PROFESSIONAL STETHOSCOPE- MIDNIGHT BLACK

602BK Deluxe Sprague-Rappaport Type Professional Stethoscope- Midnight Black
Black metal binaural
Accessory pouch included
Latex Free



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DIAGNOSTIC INSTRUMENTS](#), [Stethoscopes](#)

Tags: [professional stethoscope](#), [stethoscope](#)

1

PRODUCT DESCRIPTION

602BK Deluxe Sprague-Rappaport Type Professional Stethoscope- Midnight Black
Five interchangeable black chestpieces (2 diaphragms, 3 bells) facilitates the detection of high, medium, and low frequency sounds
Includes 3 black eartip sets to maximize comfort
Heavy-walled 22" tubing decreases extraneous sounds
Ideal for infant, pediatric, and adult use
Black metal binaural
Accessory pouch included
Latex Free

