



6433A LUMEX SILVER COLLECTION STEEL DROP ARM 3-IN-1 COMMODORE

6433A Lumex Silver Collection Steel Drop Arm 3-in-1 Commode
Arms drop independently below seat level for safe lateral transfers



SKU: N/A

Price:

Stock: 5 in stock (can be backordered)

Categories: [BATHROOM SAFETY](#), [Commodes](#)

Tags: [3-in-1 commode](#), [drop arm commode](#)

1

PRODUCT DESCRIPTION

6433A Lumex Silver Collection Steel Drop Arm 3-in-1 Commode

Arms drop independently below seat level for safe lateral transfers

Easy-to-release drop arm mechanism may be operated by patients lacking manual dexterity

Arms automatically lock in raised position

Removable backbar facilitates use as a bedside commode, raised toilet seat, or toilet safety frame

All welded steel construction increases strength and stability

Durable, grey-coated steel frame is attractive and easy to clean

Cushioned armrests enhance user comfort

Sturdy plastic seat snaps-on easily

Seat height readily adjusts from 18 $\frac{3}{4}$ " to 22 $\frac{3}{4}$ ", in 1 inch increments

Comes complete with commode pail, cover and splash shield

Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED



Designed to meet the requirements of HCPCS code: E0165
Limited lifetime warranty

1