



6438A LUMEX IMPERIAL COLLECTION 3-IN-1 STEEL DROP ARM COMMODE, 600 LB. WEIGHT CAPACITY

6438A Lumex Imperial Collection 3-in-1 Steel Drop Arm Commode
600 lb. Weight Capacity
Arms drop independently below seat level for safe lateral transfers



SKU: N/A

Price:

Stock: 3 in stock (can be backordered)

Categories: [BATHROOM SAFETY](#), [Commodes](#)

Tags: [commode](#), [drop arm commode](#)

1

PRODUCT DESCRIPTION

6438A Lumex Imperial Collection 3-in-1 Steel Drop Arm Commode
600 lb. Weight Capacity

- Arms drop independently below seat level for safe lateral transfers
- Easy-to-release mechanism may be operated by patients lacking manual dexterity
- Arms automatically lock in raised position
- Extra wide platform seat is comfortable for long periods of sitting
- Welded steel construction and extra wide base provide added stability
- Cushioned armrests enhance user comfort
- Height readily adjusts from 19 ½" - 23", in 1 inch increments
- Removable backbar facilitates use as a bedside commode raised toilet seat or toilet safety frame



- Comes complete with commode pail, cover and splash shield
- Maximum Weight Capacity 600 lb - EVENLY DISTRIBUTED
- Designed to meet the requirements of HCPCS code: E0168
- Three (3) year limited warranty