



6810A LUMEX DROP ARM VERSAMODE™ - 2 REAR LOCKING CASTERS AND 7QT PAIL

6810A Lumex Drop Arm Versamode™ - 2 Rear Locking Casters and 7qt pail
Provides transport, transfer, and commode functions all in one unit



SKU: N/A

Price:

Stock: 5 in stock (can be backordered)

Categories: [BATHROOM SAFETY](#), [Commodes](#)

Tags: [commode](#), [drop arm Versamode](#)

1

PRODUCT DESCRIPTION

6810A Lumex Drop Arm Versamode™ - 2 Rear Locking Casters and 7qt pail
Provides transport, transfer, and commode functions all in one unit

- Arms drop independently below seat level for safe lateral transfers
- Cushioned armrests enhance user comfort
- Swing-away footrests can be extended or detached
- 5" rear swivel, locking casters facilitate maneuverability in tight spaces
- Padded seat cover
- Maximum Weight Capacity 300 lb - EVENLY DISTRIBUTED
- Comes unassembled
- Comes with 7qt pail and lid
- Designed to meet the requirements of HCPCS code: E0165
- Limited Warranty

