



7016A STAINLESS STEEL DELUXE IV STAND

7016A Stainless Steel Deluxe IV Stand

Height Adjustment: 51" - 86"

Product Weight: 7.5 lb

Warranty: 2 years



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [I.V Poles / Stands](#), [Patient Room](#)

Tags: [I.V pole](#), [I.V stand](#), [steel I.v stand](#)

1

PRODUCT DESCRIPTION

7016A Stainless Steel Deluxe IV Stand

7016A IV POLE STAINLESS STEEL 4HOOK

Constructed of chrome-plated stainless steel

5 leg IV stand with 2" nylon swivel casters provides easy movement in various patient care settings

Telescoping pole adjusts from 51" - 86"

4 circular hooks for solution suspension

2 year limited warranty

Designed to meet the requirements of HCPCS code: E0776

Specifications



Height Adjustment: 51" - 86"

Product Weight: 7.5 lb

Warranty: 2 years