



7103 LUMEX PLATINUM COLLECTION 3-IN-1 STEEL COMMODE GRAY

7103 Lumex Platinum Collection 3-in-1 Steel Commode Gray

Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame



SKU: N/A

Price:

Stock: 5 in stock (can be backordered)

Categories: [BATHROOM SAFETY](#), [Commodes](#)

Tags: [3-in-1 commode](#), [steel commode](#)

1

PRODUCT DESCRIPTION

7103 Lumex Platinum Collection 3-in-1 Steel Commode Gray

Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame

Welded steel construction provides strength and stability

Increased seat depth provides added stability and comfort

Durable, white-powder coated steel frame is attractive and easy to clean

Contoured waterfall armrests are designed to enhance user comfort and provide a good gripping surface

Specifications

Width Between Arms 18 ½"

Overall Width 24 ½"

Overall Depth 20"

Seat Width 14"



Seat Depth 16 1/2"

Seat Height Range 17" - 27"

Weight Capacity 300 lb EVENLY DISTRIBUTED