



7910A-1 LUMEX SHOWER TRANSPORT CHAIR

7910A-1 Lumex Shower Transport Chair
Designed for safe transport to and from the shower



SKU: N/A

Price:

Stock: 5 in stock (can be backordered)

Categories: [BATHROOM SAFETY](#), [Shower Transport Chairs](#)

Tags: [shower chair](#), [shower transport chair](#)

1

PRODUCT DESCRIPTION

7910A-1 Lumex Shower Transport Chair

Designed for safe transport to and from the shower

Anodized, non-corrosive aluminum frame construction will not rust

Comfortable nylon backrest removes for easy cleaning

Cushioned armrests enhance user comfort

5" swivel casters with rear wheel safety locks facilitate maneuverability in tight spaces

Rolls easily over most toilets

Allows complete access to body for easy cleaning during bathing

Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED

Limited lifetime warranty

1