



## 7915A-1 DELUXE SHOWER TRANSPORT CHAIR WITH FOOTRESTS

7915A-1 Deluxe Shower Transport Chair with Footrests  
Freight Exception Aluminum  
Designed for safe transport to and from the shower



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [BATHROOM SAFETY](#), [Shower Transport Chairs](#)

**Tags:** [shower transport chair](#), [transport chair](#)

1

### PRODUCT DESCRIPTION

7915A-1 Deluxe Shower Transport Chair with Footrests

Freight Exception Aluminum

Designed for safe transport to and from the shower

Anodized, non-corrosive aluminum frame construction will not rust

Comfortable nylon backrest removes for easy cleaning

Cushioned armrests enhance user comfort

5" swivel casters with rear wheel safety locks facilitate maneuverability in tight spaces

Rolls easily over most toilets

Allows complete access to body for easy cleaning during bathing

Powder-coated steel front brace with footrest

Foot plates tip up and out of the way for unobstructed entry and exit

Footrests comfortably support the feet during transport



Maximum Weight Capacity 250 lb - EVENLY DISTRIBUTED  
Limited lifetime warranty