



80072 CLINTON IMAGING TABLE W FOWLER BACK AND DROP WINDOW

80072 Clinton Imaging Table W Fowler Back And Drop Window
Item No: 80072



SKU: N/A

Price:

Stock: Available on backorder

Categories: [IMAGING TABLES](#), [Ultrasound Table](#)

Tags: [imaging table](#), [ultrasound table](#)

1

PRODUCT DESCRIPTION

80072 Clinton Imaging Table W Fowler Back And Drop Window
Item No: 80072

CLINTON TWO-PIECE TOP WITH PNEUMATIC ADJUSTABLE BACKREST AND WINDOW DROP; TRENDLENBURG/REVERSE TRENDLENBURG PLUS OR MINUS 15; EASY-TO-USE HAND CONTROL; ALL-WELDED BASE AND TOP FRAME; UNOBSTRUCTED BASE DESIGN; DUAL COLUMN LIFT SYSTEM ADDS RIGIDITY; PAPER DISPENSER INCLUDED; POWER SUPPLIES AND MOTORS MEET UL 60601-1 AND CAN/CSA C22.2 601.1; INDUSTRY LEADING 5-YEAR, LIMITED WARRANTY; 34 WIDE MODEL AVAILABLE; 550 LBS LOAD CAPACITY UNDER NORMAL USE

1