



## 80073 CLINTON IMAGING TABLE WITH 3-SECTION TOP AND DROP WINDOW

80073 Clinton Imaging Table With 3-Section Top And Drop Window

Item No: 80073



**SKU:** N/A

**Price:**

**Stock:** Available on backorder

**Categories:** [IMAGING TABLES](#), [Ultrasound Table](#)

**Tags:** [imaging table](#), [ultrasound table](#)

1

### PRODUCT DESCRIPTION

80073 Clinton Imaging Table With 3-Section Top And Drop Window

Item No: 80073

CLINTON THREE-PIECE TOP WITH PNEUMATIC ADJUSTABLE BACKREST, MANUAL FOOTREST AND WINDOW DROP; FOOTREST IS ADJUSTABLE FROM EITHER SIDE AND RAISES UP TO 30 FROM FLAT POSITION; OPTIONAL SIDE RAILS & CASTERS (SEE OPT); TRENDLENBURG/REVERSE TRENDLENBURG PLUS OR MINUS 15; EASY-TO-USE HAND CONTROL; ALL-WELDED BASE AND TOP FRAME; UNOBSTRUCTED BASE DESIGN; DUAL COLUMN LIFT SYSTEM ADDS RIGIDITY; PAPER DISPENSER INCLUDED; POWER SUPPLIES AND MOTORS MEET UL 60601-1 AND CAN/CSA C22.2 601.1; INDUSTRY LEADING 5-YEAR, LIMITED WARRANTY; 34 WIDE MODEL AVAILABLE;

1