



85100 CLINTON POWER IMAGING TABLE W/ WINDOW DROP OPEN BASE

85100 Clinton Power Imaging Table W/ Window Drop Open Base

Made in the USA by Clinton Industries

450 lbs. load capacity under normal use



SKU: N/A

Price:

Stock: Available on backorder

Categories: [IMAGING TABLES](#), [Ultrasound Table](#)

Tags: [imaging table](#), [ultrasound table](#)

1

PRODUCT DESCRIPTION

85100 Clinton Power Imaging Table W/ Window Drop Open Base

Item No: 85100

Features flat top with 10" x 10" hinged imaging window section with metal slide latch

Electric power height adjustment with foot control

Ultra-strong, uni-frame base and top construction

Heavy duty, all welded, lift mechanism

Durable gray powder-coated steel frame

Self-lubricating pivot points

Smooth 110 volt electric height adjustment

Power supply and motor meets UL 60601-1 and CAN/CSA C22.2 601.1

Hands-free low voltage foot control

Easy service, self contained, electrical components



Rounded corner, upholstered top
2" ultra-firm padding
Low height may help achieve ADA tax credit
Customize with options
Made in the USA by Clinton Industries
450 lbs. load capacity under normal use

1