



89100-KD LUMEX 18" PVC KNOCK-DOWN SHOWER COMMODE CHAIR

89100-KD Lumex 18" PVC Knock-Down Shower Commode Chair

Easy to assemble and breakdown for convenient transport

Extensively tested



SKU: N/A

Price:

Stock: 4 in stock (can be backordered)

Categories: [BATHROOM SAFETY](#), [Shower Transport Chairs](#)

Tags: [shower chair](#), [shower commode](#)

STANDARD MESH COLORS



R. Blue



F. Green



Mauve

1

SOLID VINYL COLORS



R. Blue



F. Green



Mauve

GALLERY IMAGES





PRODUCT DESCRIPTION

89100-KD Lumex 18" PVC Knock-Down Shower Commode Chair

- Easy to assemble and breakdown for convenient transport
- Constructed of the highest quality health care grade PVC
- Deluxe elongated extra heavy-duty, contoured open front seat
- Dual reinforcement into the caster fittings adds additional support
- Quick drying, removable mesh back support is standard on all shower chairs
- Available in 3 anti-bacterial mesh or anti-bacterial solid vinyl colors
- Ergonomic curved arms provide support without sharp edges and include anti-slip handgrips for added safety
- Threaded stem casters for additional safety
- One year limited warranty on PVC frame; six month limited warranty on all other materials
- Clamp Joint Bracket Locking System is made of high grade impact plastic and stainless steel
- No sharp edges
- Will secure pipe and fittings 100%
- No tools required
- Extensively tested
- Patent pending on clamp joint bracket