



1307 BABINSKI PERCUSSION HAMMER

1307 Babinski Percussion Hammer
1307 9" with chrome-plated handle.



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [babinski hammer](#), [percussion hammer](#)

1

PRODUCT DESCRIPTION

1307 Babinski Percussion Hammer
1307 9" with chrome-plated handle.
Adjustable head can be oriented in two directions
Hammer contains neuro monofilament for neurologic testing

1