



## BARIATRIC TRAPEZE

Bariatric Trapeze

2940B Lumex Bariatric Trapeze, 450 lb weight capacity

2960B Lumex Bariatric Trapeze, 600 lb weight capacity



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [Patient Room](#), [Trapez / Floor Stand](#)

**Tags:** [bariatric trapeze](#), [patient trapeze](#)

### ADDITIONAL INFORMATION

**Item No. 2**                      281, 282

### VARIATIONS

**Image**

**SKU**

**Price**


**Description**

**Item No. 2**



2940B



Image	SKU	Price	Description	Item No. 2
				2960B

**PRODUCT DESCRIPTION**

**Bariatric Trapeze**

2940B Lumex Bariatric Trapeze, 450 lb weight capacity

2960B Lumex Bariatric Trapeze, 600 lb weight capacity

Helps patient safely change position in bed and transfer from bed to other equipment with minimal assistance.

Rubber bumper guards protect walls

Handbar location fully adjustable

Comes with floor stand

Meets requirements for HCPCS code E0912

14mm wrench required for assembly

**Specifications 2940B 2960B**

Assisted Weight Capacity 450 lb 600 lb

Overall Height (including fixed bracket) 76 7/8" 80 5/8"

Outside Base Width 28 7/8" 29 1/4"

Width Between Base Legs 27" 25 1/4"

Base Length 40 1/2" 41"

Weight of Product 72 lbs 90 lbs

**ADDITIONAL INFORMATION**

<b>Item No. 2</b>	281, 282
-------------------	----------

