



## BATTERY-POWERED FLOOR

Battery-Powered Floor Lifts

LF1050 400 lb weight capacity

DPL650BAT BATTERY LF1050 1090 2020 2090

DPL650PEN HAND CTL LF1050 1090 2020 2090

DPL650ADP CHARGER LF1050 1090 2020 2090



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [Battery Powered Lifts](#), [PATIENT LIFTS](#)

**Tags:** [battery-powered lift](#), [floor lift](#), [patient lift](#)

## ADDITIONAL INFORMATION

**Item No. 2**                      367, 365, 366, 370

## VARIATIONS



**Image**

**SKU**

**Price**

**Description**

**Item No. 2**



LF1050

**PRODUCT DESCRIPTION**

**Battery-Powered Floor Lifts**

LF1050 400 lb weight capacity

DPL650BAT BATTERY LF1050 1090 2020 2090

DPL650PEN HAND CTL LF1050 1090 2020 2090

DPL650ADP DPL650ADP-CABLE Replacement Adapter Cable for Charger LF1050 1090 2020 2090

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Six point spreader bar with 360° rotation accommodates 2-, 4-, or 6-point slings
- Large, easy to grip handles for easy maneuverability
- Intuitive foot pedal leg opening mechanism frees up hands
- Ergonomic hand pendant with hook
- Emergency stop button and manual lowering mechanism
- Warning beep and visual indicator on control box for low battery charge (24V DC motor)
- Accessory storage bag and how-to-use video included
- Optional digital scale DSC250
- Designed to meet the requirements of HCPCS Code E0635

**Warranty:**

Frame and spreader bar: Three (3) year limited

Actuator: Two (2) year limited

Batteries: One (1) year limited



**ADDITIONAL INFORMATION**

**Item No. 2** 367, 365, 366, 370