



BATTERY-POWERED SIT-TO-STAND LIFTS

Battery-Powered Sit-to-Stand Lifts

LF2020 LIFT SIT TO STAND 400 LB

DPL650BAT BATTERY LF1050 1090 2020 2090

DPL650PEN HAND CTL LF1050 1090 2020 2090

DPL650ADP CHARGER LF1050 1090 2020 2090



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Battery Powered Lifts](#), [PATIENT LIFTS](#)

Tags: [battery-powered lift](#), [patient lift](#), [sit-to-stand lift](#)

ADDITIONAL INFORMATION

Item No. 2 367, 365, 366, 371

VARIATIONS



Image	SKU	Price	Description	Item No. 2
				LF2020

PRODUCT DESCRIPTION

Battery-Powered Sit-to-Stand Lifts

LF2020 LIFT SIT TO STAND 400 LB

DPL650BAT BATTERY LF1050 1090 2020 2090

DPL650PEN HAND CTL LF1050 1090 2020 2090

DPL650ADP CHARGER LF1050 1090 2020 2090

Sit-to-Stand lifts are designed for patients who have difficulty standing up on their own from a seated position.

Patients must be able to actively participate in the lifting process and be able to partially support their weight.

- Heavy gauge steel construction with white powder-coated finish and rubber-coated low base legs protect furniture and walls
- Spreader bar provides various hand-hold locations for the patient and accommodates various sling styles
- Large, easy to grip handles for easy maneuverability
- Intuitive foot pedal leg opening mechanism frees up hands
- Adjustable knee pad unit
- Quick-release removable foot platform for gait training
- Ergonomic hand pendant with hook
- Emergency stop button
- Warning beep and visual indicator on control box for low battery charge (24V DC motor)
- Accessory storage bag and how-to-use video included
- Designed to meet the requirements of HCPCS Code E0635



Warranty:

Frame and spreader bar: Three (3) year limited

Actuator: Two (2) year limited

Batteries: One (1) year limited

ADDITIONAL INFORMATION

Item No. 2 367, 365, 366, 371