



BED SOCKET TELESCOPING I.V. POLE

Bed Socket Telescoping I.V. Pole
 2100A IV Pole Bed Socket Telescoping 3 / CS
 2101A IV Pole Bed Socket Telescoping INDIVIDUALL



SKU: N/A
Price:
Stock: 10 in stock (can be backordered)
Categories: [I.V Poles / Stands](#), [Patient Room](#)
Tags: [I.V pole](#), [medical stand](#), [telescoping I.V pole](#)

ADDITIONAL INFORMATION

Package 1615, 1692

VARIATIONS



Image	SKU	Price	Description	Package
				Case



Image	SKU	Price	Description	Package
				Individuall

PRODUCT DESCRIPTION

Bed Socket Telescoping I.V. Pole

2100A IV Pole Bed Socket Telescoping 3 / CS

2101A IV Pole Bed Socket Telescoping INDIVIDUALL

- Locking collet quickly and securely adjusts height from 35" to 65".
- Internal rubber bumpers ensure silent operation and reduce shock to hand and I.V. container when lowering telescoping tube.
- Constructed of durable, anodized aluminum tubing.
- Designed to meet the requirements of HCPCS code: E0776.

ADDITIONAL INFORMATION

Package 1615, 1692