



CARDIOLOGY STETHOSCOPE

Cardiology Stethoscope
425 Cardiology Stethoscope Grey
425BLK Cardiology Stethoscope Black



SKU: N/A

Price:

Stock: 18 in stock

Categories: [DIAGNOSTIC INSTRUMENTS](#), [Stethoscopes](#)

Tags: [cardiology stethoscope](#), [stethoscope](#)

ADDITIONAL INFORMATION

Weight	1.00 lbs
Dimensions	14 × 8 × 3 in
Color	107, 110

VARIATIONS



Image	SKU	Price	Description	Color
				Black
				Grey

PRODUCT DESCRIPTION

Cardiology Stethoscope
 425 Cardiology Stethoscope Grey
 425BLK Cardiology Stethoscope Black
 Features a stainless steel combination chestpiece for maximum convenience and internal spring binaurals for high acoustic sensitivity. Dual non-chill rings provide patient comfort and 22" latex-free tubing keeps acoustics clear, even in the noisiest settings. Comes complete with additional diaphragm and two pair of contoured eartips. Lifetime Warranty.

ADDITIONAL INFORMATION

Weight	1.00 lbs
Dimensions	14 × 8 × 3 in
Color	107, 110

