




## CHROME PLATED EAR SYRINGE

Chrome Plated Ear Syringe  
 2807 2oz. Ear Syringe  
 2808 3oz. Ear Syringe  
 2809 4oz. Ear Syringe  
 GRA2810 Shield Only



**SKU:** N/A

**Price:**

**Stock:** 10 in stock

**Categories:** [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

**Tags:** [chrome ear syringe](#), [plated ear syringe](#)

### ADDITIONAL INFORMATION

**volume** 117, 118, 119, 121

### VARIATIONS





Image	SKU	Price	Description	volume
				2oz



Image	SKU	Price	Description	volume
				3oz
				4oz
				Shield Only

## PRODUCT DESCRIPTION

Chrome Plated Ear Syringe

2807 2oz. Ear Syringe

2808 3oz. Ear Syringe

2809 4oz. Ear Syringe

GRA2810 Shield Only

Grafco Ear Syringes feature a durable brass construction with attractive chrome plating. Ear Irrigation Syringes are ideal for washing and gently clearing ear wax and debris. The Ear Syringes are completely latex free. Available in 3 different sizes, each Ear Syringe includes a shield and two tips.

- Brass construction
- Chrome plated
- Includes shield and two tips
- Available in 3 sizes
- Latex free
- 1 Year Limited Warranty: Product is warranted to be free from defects in material and workmanship. Repair or replacement will be issued for any product which is found defective under the terms of this warranty.



Specifications

Item#	Unit I.D	Unit Length	Tip I.D.	Tip Length
2807	1.122"	6.625"	0.0890" & 0.0940"	1.906" & 1.998"
2808	1.403"	6.750"	0.0940" & 0.0960"	1.877" & 1.951"
2809	1.403"	7.750"	0.0940" & 0.0975"	1.870" & 2.003"

**ADDITIONAL INFORMATION**

**volume** 117, 118, 119, 121