




COMPLETE DIAGNOSTIC SET

Complete Diagnostic Set
1227-4 Complete Diagnostic Set
1229-1 Bayonet-Catch Adapter For Disposable Specula
1230L Ophthalmoscope Bulb
1230X 9 mm Aural/Nasal Speculum
1230Y 1.25 mm Aural Speculum



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL and PHYSICIAN](#)

Tag: [diagnostic set](#)

1

PRODUCT DESCRIPTION

Complete Diagnostic Set
1227-4 Complete Diagnostic Set
1229-1 Bayonet-Catch Adapter For Disposable Specula (required for use with Item 1229 when using Items 1224, 1226, 1227-4)
1230L Ophthalmoscope Bulb (for Items 1227-2, 1227-4, 355)
1230X 9 mm Aural/Nasal Speculum (For use with Items 1224, 1226, 1227-1, 1227-4, 1228-3, 5100X)
1230Y 1.25 mm Aural Speculum (For use with Items 1224, 1226, 1227-1, 1227-4, 1228-3,



5100X)

Complete set for eyes, ears, nose and throat.

Includes:

Chrome-plated otoscope head with 3 aural bayonet-catch specula (small, medium, large)

Ophthalmoscope head with metal nasal speculum

Metal tongue depressor and throat lamp

#2 18mm and #4 22mm mirror attachments

Chrome-plated screw lock handle with rheostat houses

Replacement bulb and case

Sterilizable and autoclavable. Houses 2 C Batteries (Not Included).