



DUAL HEAD STETHOSCOPE

Dual Head Stethoscope

400 Black

400BL Blue

400GR Green

400GY Grey

400L Lavender

400P Hot Pink

400R Red



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [DIAGNOSTIC INSTRUMENTS](#), [Stethoscopes](#)

Tags: [dual head stethoscope](#), [stethoscope](#)

ADDITIONAL INFORMATION

Color 107, 251, 256, 110, 250, 257, 85

VARIATIONS



Image

SKU

Price

Description

Color



Red



Black



Grey



Lavender



Image **SKU** **Price** **Description** **Color**



Blue



Green



Pink

PRODUCT DESCRIPTION

Dual Head Stethoscope

400 STETH DUAL HEAD BLK Black

400BL STETH DUAL HEAD BL Blue

400GR STETH DUAL HEAD GREEN Green

400GY STETH DUAL HEAD GREY Grey

400L STETH DUAL HEAD LAVENDER Lavender

400P STETH DUAL HEAD PINK Hot Pink

400R STETH DUAL HEAD RED Red

Anodized aluminum rotating chestpiece with recessed non-chill ring and adjustable steel



binaural.

22" 'Y' tubing not made with natural rubber latex

ADDITIONAL INFORMATION

Color 107, 251, 256, 110, 250, 257, 85