



EMESIS BASINS

Emesis Basins

3231 Small Capacity: 13 oz.

3232 Large Capacity: 26 oz.



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [MEDICAL, SURGICAL and PHYSICIAN, Stainless Steelware](#)

Tags: [emesis basins](#), [medical basins](#)



ADDITIONAL INFORMATION

Size 90, 88

VARIATIONS

Image	SKU	Price	Description	Size
-------	-----	-------	-------------	------



Image	SKU	Price	Description	Size
				S
				L

PRODUCT DESCRIPTION

Emesis Basins

3231 Small, 8 1/8" x 3 5/8" x 1 3/4". Capacity: 13 oz.

3232 Large, 10" x 4 1/2" x 2 1/8". Capacity: 26 oz.

18-8 (304) Stainless Steel with seamless construction.

ADDITIONAL INFORMATION

Size 90, 88