



FAST ALERT ADVANCED PATIENT ALARM

Fast Alert Advanced Patient Alarm


GF13702 Lumex Fast Alert Advanced Patient Alarm

GF13702B Lumex Fast Alert Advanced Patient Alarm with Bed Pad

GF13702C Fast Alert Advanced Patient Alarm with Chair Pad

106375 Pressure Sensor Bed Pad (1 Year Warranty)

106380 Pressure Sensor Chair Pad (1 Year Warranty)



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Patient Alarms](#), [Patient Room](#)

Tags: [advanced alarm](#), [alert alarm](#), [patient alarm](#)

ADDITIONAL INFORMATION

Item No. 2 332, 333, 329, 330, 331

VARIATIONS



Image	SKU	Price	Description	Item No. 2
				GF13702
				GF13702B
				GF13702C
				106375



Image

SKU

Price

Description

Item No. 2



106380

PRODUCT DESCRIPTION

Fast Alert Advanced Patient Alarm

GF13702 Lumex Fast Alert Advanced Patient Alarm White

GF13702B Lumex Fast Alert Advanced Patient Alarm with Bed Pad White

GF13702C Fast Alert Advanced Patient Alarm with Chair Pad White

106375 Pressure Sensor Bed Pad (1 Year Warranty) White

106380 Pressure Sensor Chair Pad (1 Year Warranty) White

Alarm notifies staff or caregiver when user gets out of bed or chair

Pressure sensitive pad with 1 year warranty connects to alarm

Alarm easily secures with self contained clip onto areas such as wheelchairs

Convenient dual volume settings: 70dB and 90dB

Delay function allows for 0 or 2 second delay to reduce likelihood of false alarms with restless users

Reset button on side of monitor silences and resets alarm

Remote alarm option allows the caregiver to silence the local alarm and sound the alarm remotely through the nurse call system (Adapter Cable and Splitter sold separately)

Requires one 9V battery (not included)

Comes with protective boot cover to prevent damage if dropped

Off switch and Low battery indicator

1 Year Limited Warranty

Pad size: Chair - 10" x 15"

Pad size: Bed - 10" x 29"



ADDITIONAL INFORMATION

Item No. 2 332, 333, 329, 330, 331