



FAST ALERT BASIC PATIENT ALARM

Fast Alert Basic Patient Alarm

GF13701 Lumex Fast Alert Basic Patient Alarm

GF13701B Lumex Fast Alert Basic Patient Alarm with Bed Pad

GF13701C Lumex Fast Alert Basic Patient Alarm with Chair Pad

GF13701CR PATIENT ALARM BASIC W/CHAIR PAD RETAIL PACKAGED

GF13701BR PATIENT ALARM BASIC W/ BED PAD RETAIL PACKAGED

106375 Pressure Sensor Bed Pad (1 Year Warranty)

106380 Pressure Sensor Chair Pad (1 Year Warranty)



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Patient Alarms](#), [Patient Room](#)

Tags: [alert alarm](#), [patient alarm](#)

ADDITIONAL INFORMATION

Item No. 2 332, 333, 338, 339, 342, 340, 341

VARIATIONS



Image

SKU Price Description

Item No. 2



GF13701



GF13701B



GF13701C



GF13701CR



Image **SKU** **Price** **Description** **Item No. 2**



GF13701BR



106375



106380

PRODUCT DESCRIPTION

Fast Alert Basic Patient Alarm

GF13701 Lumex Fast Alert Basic Patient Alarm White

GF13701B Lumex Fast Alert Basic Patient Alarm with Bed Pad White

GF13701C Lumex Fast Alert Basic Patient Alarm with Chair Pad White

GF13701CR PATIENT ALARM BASIC W/CHAIR PAD RETAIL PACKAGED

GF13701BR PATIENT ALARM BASIC W/ BED PAD RETAIL PACKAGED

106375 Pressure Sensor Bed Pad (1 Year Warranty) White

106380 Pressure Sensor Chair Pad (1 Year Warranty) White

Alarm notifies staff or caregiver when user gets out of bed or chair



Pressure sensitive pad with 1 year warranty connects to alarm
Alarm easily secures with self contained clip onto areas such as wheelchairs
Convenient dual volume settings: 70dB and 90dB
Delay function allows for 0 or 2 second delay to reduce likelihood of false alarms with restless users
Requires one 9V battery (not included)
Comes with protective boot cover to prevent damage if dropped
Off switch and Low battery indicator
1 Year Limited Warranty
Pad size: Chair - 10" x 15"
Pad size: Bed - 10" x 29"
GF1370-COVER Replacement Protective Boot Cover for Patient Alarm
GF1370-BCLIP Bed Clip for Patient Alarm

ADDITIONAL INFORMATION

Item No. 2 332, 333, 338, 339, 342, 340, 341