




GF7012-1 SELECT CARE 2-HOOK I.V. STAND

GF7012-1 Select Care 2-Hook I.V. Stand

Base width: 27.5" x 27.5"

Height adjustable from 47" - 84"

Designed to meet the requirements of HCPCS code: E0776



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [I.V Poles / Stands](#), [Patient Room](#)

Tags: [2-hook stand](#), [I.V stand](#), [medical stand](#), [stands](#)

1

PRODUCT DESCRIPTION

GF7012-1 Select Care 2-Hook I.V. Stand

Low profile base prevents tipping

Swivel ball bearing casters facilitate a smooth transport of the pole on all surfaces

Base width: 27.5" x 27.5"

Height adjustable from 47" - 84"

Chrome plated steel material

Designed to meet the requirements of HCPCS code: E0776

1