



GF8902 ECONOMY OVERBED TABLE, NON-TILT

GF8902 Economy Overbed Table, Non-Tilt

GF8902 Lumex Everyday Overbed Table, Non-Tilt Wood & Chrome plated

Easy to assemble



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Overbed Table](#), [Patient Room](#)

Tags: [non-tilt table](#), [overbed table](#)

1

PRODUCT DESCRIPTION

GF8902 Economy Overbed Table, Non-Tilt

GF8902 Lumex Everyday Overbed Table, Non-Tilt Wood & Chrome plated

Easy to assemble

Height adjusts infinitely between 28" and 44"

Twin-wheel swivel casters provide smooth-rolling action. Two locking and two non-locking casters

Walnut, wood-grain laminated 3/4" particle-board core with plastic edge molding

Chrome plated "H" base and column for aesthetics and durability

Designed to meet the requirements of HCPCS code: E0274

Supports up to 30 lbs. of dispersed weight.

Warranty : One Year

Specifications



Table Height 28 - 44"

Table Width 30"

Table Depth 15"

Base Width 26.2"

Base Depth 15.25"

Weight 22 lbs

Dispersed Weight Capacity 30 lbs