



## GF8905-1A DELUXE TILT OVERBED TABLE

GF8905-1A Deluxe Tilt Overbed Table

Designed to meet the requirements of HCPCS code: E0274.

Supports up to 30 lb of dispersed weight.

One year warranty.



**SKU:** N/A

**Price:**

**Stock:** 10 in stock (can be backordered)

**Categories:** [Overbed Table](#), [Patient Room](#)

**Tags:** [overbed table](#), [tilt overbed](#)

1

## PRODUCT DESCRIPTION

GF8905-1A Deluxe Tilt Overbed Table

Easy assembly.

Chrome-plated, welded tubular steel column and base provide durability and stability.

1 ½" hooded, ball-bearing casters ensure easy maneuverability across all floor surfaces.

Table top, constructed of wood-grained oak Melamine laminate on particle-board core, adjusts 13" in infinite positions.

Raised T-molded plastic edge.

Chrome-plated "H" base and column for aesthetics and durability.

2 locking casters.

Designed to meet the requirements of HCPCS code: E0274.

Supports up to 30 lb of dispersed weight.



One year warranty.

1