



GP 300 INFRARED MEDICAL NON-CONTACT THERMOMETER

GP 300 Infrared Medical Non-Contact Thermometer


Product Details:

Infrared Technology

Contactless

One second temperature measurement.

Memory: Memorize 32 measured values



SKU: N/A

Price:


Stock: 15 in stock (can be backordered)

Category: [PPE](#)

1



GALLERY IMAGES





PRODUCT DESCRIPTION

GP 300 Infrared Medical Non-Contact Thermometer

Product Description

It is the best product for taking the temperature of the body and surfaces.

Backed by a team of well-informed experts, we are a well-known firm in the market for providing an extensive series of GP 300 Taiols Non-Contact Infrared [Thermometer](#). Moreover, we offer this product at nominal costs.

This is a contactless forehead temperature with a measurement range of 32C-42.9C (89.6F-109.2F), measurement accuracy within (35-42)C. Automatic shutdown without any operation within 30 seconds, product size is 95mm x 60mm x 180mm. WEIGHT OF PRODUCT IS 146g. Memory

Product Details:

Infrared Technology

Contactless

One second temperature measurement.

Memory: Memorize 32 measured values