



GRAFCO® LARYNGOSCOPE SETS


Grafco Laryngoscope Sets

1196 MacIntosh Standard Laryngoscope Set

1196F MacIntosh Fiberoptic Laryngoscope Set

1197 Miller Standard Laryngoscope

1197F Miller Fiberoptic Laryngoscope Set



SKU: N/A

Price:

Stock: 10 in stock (can be backordered)

Categories: [Instruments](#), [MEDICAL](#), [SURGICAL](#) and [PHYSICIAN](#)

Tags: [grafco sets](#), [laryngoscope sets](#)

ADDITIONAL INFORMATION

Item No. 2 128, 127, 130, 129

VARIATIONS



Image	SKU	Price	Description	Item No. 2
				1196f
				1196
				1197f
				1197

PRODUCT DESCRIPTION

Grafco Laryngoscope Sets

1196 MacIntosh Standard Laryngoscope Set: Includes 1 Medium "C" battery handle, 4 blades



(1 ea. #1, #2, #3, #4) and carry case

1196F MacIntosh Fiberoptic Laryngoscope Set: Includes 1 Medium "C" battery handle, 4 blades (1 ea. #1, #2, #3, #4) and carry case

1197 Miller Standard Laryngoscope Set: Includes 1 Medium "C" battery handle, 4 blades (1 ea. #1, #2, #3, #4) and carry case

1197F Miller Fiberoptic Laryngoscope Set: Includes 1 Medium "C" battery handle, 4 blades (1 ea. #1, #2, #3, #4) and carry case

Stainless steel construction with non-glare satin finish. Set includes 1 Medium "C" battery handle, 4 blades (1 ea. #1,2,3,4) and handy carry case.

ADDITIONAL INFORMATION

Item No. 2 128, 127, 130, 129